


# Pathfinder 164

#### **Parameters**

Model	Usable Length	Profile	Lens Position From Tip	Separate Lens Lumen
Pathfinder 164	1350mm	2. 58Fr(0.86mm)	20mm	Yes

#### Excellent lesion crossing ability



Low-profile: 2. 58Fr(0.86mm)



Hydrophobic coating



Enhanced out-layer of distal catheter

### Outstanding capability to observe more distant artery



Short distance from Lens to Tip: 20mm



## Shenzhen Vivolight Medical Device & Technology Co., Ltd.

Address: 2F-5F, GAOXINQIKEJILOU Building, Block 67, Xingdong community, Xin'an street,

Bao'an District, Shenzhen, P.R. China

Email: sales@vivo-light.com

Website: www.vivolight.com