

Pathfinder 164

Parameters

Model	Usable Length	Profile	Lens Position From Tip	Separate Lens Lumen
Pathfinder 164	1350mm	2. 58Fr(0.86mm)	20mm	Yes

Excellent lesion crossing ability



Low-profile: 2. 58Fr(0.86mm)



Hydrophobic coating



Enhanced out-layer of distal catheter

Outstanding capability to observe more distant artery



Short distance from Lens to Tip: 20mm



Shenzhen Vivolight Medical Device & Technology Co., Ltd.

Address: 2F-5F, GAOXINQIKEJILOU Building, Block 67, Xingdong community, Xin'an street,

Bao'an District, Shenzhen, P.R. China

Email: sales@vivo-light.com **Website:** www.vivolight.com